

CellSonic Ltd

www.cellsonic-medical.com

www.cellsonic-beauty.com





Synchrowave

Synchrowave is a physiotherapy machine which synchronises the action of a radial shockwave with an electrotherapy current to give enhanced muscle stimulation. At the same time that the radial shockwave fires the electrotherapy pulses and the muscle is both tensed and massaged. The result is quicker improvement in body function and better pain control.

Position the electro-tabs either side of the area to be treated and aim the Synchrowave hand-piece between the tabs keeping it moving slowly from pulse to pulse.

The Synchrowave runs off compressed air and the compressor supplied by CellSonic is noted for being very quiet.



**Quality approved to ISO
13485:2003 & ISO 9001:2008**

CE 1023

All enquiries to sales@cellsonic-medical.com